

SMARTFIT

CAD/CAM CUSTOM ABUTMENT



SMARTFIT is Hiossen's solution to your custom abutment needs. Hiossen utilizes the most state of the art Computer Aided Design/Computer Aided Manufacturing (CAD/CAM) technology available to fabricate a truly unique, one of a kind abutment for your patient.

Hiossen's SMARTFIT team has masterfully combined art and technology to create the perfect abutment, for the perfect smile.



*Zirconia & Bar Type custom abutments will be available in 2015


OSSTEM[®]
IMPLANT

PROCESS OVERVIEW

1. Preparation



Implant placement



Fixture level impression



Working model

2. Scan & Design



Digital Scan of working model



3D Scanned image



SMARTFIT custom abutment design

3. Milling & Fitting



Milled SMARTFIT custom abutment



SMARTFIT w/ a jig



Crown fabrication

4. Placement & Final Restoration



SMARTFIT connected w/ jig



SMARTFIT Final



Final crown connected

To order a SMARTFIT custom abutment, please contact your local sales representative or the SMARTfit department: